ANYCHECK

*2025 Global Top Dental Implant Company


IF you want to check Loading time...

ANYCHECK

Anycheck :

- 何謂 Loading?

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
The above is a summary of scientific data and not an official recommendation by Osstell.
To monitor osseointegration measure at placement and before final restoration.

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(Anycheck vs Periotest M)

Article

Clinical Validation of Dental Implant Stability by Newly Designed Damping Capacity Assessment Device during the Healing Period

Ho-Kyung Lim ^{1,†}, Sung-Jae Lee ^{2,†}, Yujin Jeong ³, Jong-Seok Lee ⁴, Jae-Jun Ryu ⁴, Ji-Suk Shim^{5,*} and In-Seok Song ^{2,*}

¹ Department of Oral & Maxillofacial Surgery, Korea University Guro Hospital, Seoul 08308, Korea

² Department of Oral & Maxillofacial Surgery, Korea University Anam Hospital, Seoul 02841, Korea

³ Department of Biostatistics, Korea University College of Medicine, Seoul 02841, Korea

⁴ Department of Prosthodontics, Korea University Anam Hospital, Seoul 02841, Korea

⁵ Department of Prosthodontics, Korea University Guro Hospital, Seoul 08308, Korea

* Correspondence: shoss@hanmail.net (J.-S.S.); densis@korea.ac.kr (I.-S.S.)

† These authors contributed equally to this work.

Background and Objectives:

To evaluate the stability of a dental implant and the effectiveness of a newly designed damping capacity assessment device by improving the number of blows and strength evaluated by a prospective clinical study.

Materials and Method:

The stability of dental implants was measured in 50 implants in a total of 38 patients.

Measurements were performed using Anycheck and Periotest M devices, twice in total, divided into buccal and lingual directions.

In addition, measurements were performed on the day of surgery, two weeks, one month, two months, and three months after surgery for a total of five times.

After the standardization of the measured values, the differences and changes over time for each device were observed.

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Materials and Method:



Figure 1. Clinical application of Anycheck and Periostest M equipment. (A) Anycheck, (B) Periostest M).

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Result:

No difference in standardized values between the two devices was observed at any time point.

In both devices, stability decreased at two weeks postoperatively but gradually increased thereafter.

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
Figure 4. The predicted value over time (Z-score standardized, (A) buccal, (B) lingual, (C) mean).

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Conclusions:

The damping capacity of Anycheck was similar to that of Periotest M.

After a slight decrease in stability two weeks after implant placement, implant stability increased over time.

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RESEARCH

Open Access

Validation of an implant stability measurement device using the percussion response: a clinical research study

Yurie Okuhama¹, Koudai Nagata¹, Hyunjin Kim¹, Hayato Tsuruoka¹, Mihoko Atsumi² and Hiromasa Kawana^{1*}

Background and Objectives:

This study aimed to determine the correlation between primary and secondary implant stability using the **Osstell®** and **AnyCheck®** devices.

Materials and Method:


Ten patients (7 women) with partially edentulous jaws who received a total of 15 implants were included. IST (AnyCheck®) and ISQ (Osstell®) values were measured immediately after implantation and at 1, 2, 3, 4, and 6 weeks post-implantation.

Each measurement was performed three times, and the average value was used as the result.

The correlation between measurements obtained using the two devices was determined using Spearman's rank correlation coefficient.

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
Materials and Method:

Results:

The IST values ranged from 79.1 ± 2.87 to 82.4 ± 2.65 .

The ISQ values ranged from 76.0 ± 2.8 to 80.2 ± 2.35 .

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
Results:

Spearman's rank correlation coefficient was $r = 0.64$ immediately after implantation, $r = 0.29$ at 1 week, $r = 0.68$ at 2 weeks, $r = 0.53$ at 3 weeks, $r = 0.68$ at 4 weeks, and $r = 0.56$ at 6 weeks.

A positive correlation was found in all cases, except at week 1 when the correlation was weak; the IST and ISQ values decreased the most during the first postoperative week and increased during the second week.

The IST values were also slightly higher at all measurement points.

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Conclusion:

The ability to assess implant stability **without removing the abutment during healing** is essential for determining the timing of loading without the risk of bone resorption.

The results of this study suggest that **AnyCheck® is useful for determining primary and secondary implant stability.**

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