

NEOBITECH DIGITAL SYSTEM **Work Flow**

01

SCAN



Neobiotech
Scan Body
IS/IT/MUA用



02

CAD/CAM



Neobiotech Digital
Libraries Available
on Exocad, 3shape,
and Dental wings

03

Milling



Neobiotech
Premill Ti-Block
IS用 鈦棒



04 Custom abutment Materials

IS用 Ti-Base Components

Anterior : Ti-Link



Posterior : Ti-Cylinder



IT用 Ti-Link




MUA用 Cylinder



05 Furnace



06 Surgical guide



IS用 Lab Analog




IT用 Lab Analog



Surgical Guide Sleeve



Scan Body



IS / IT用




Implant	Type	Ø4.0	Ø5.0
Neobiotech IS	Hex	ISPSBH40NB	ISPSBH50NB
	SCRIP	ISPSBS40NB	ISPSBS50NB



Implant	Type	Ø4.8	Ø6.5
Neobiotech IT	Octa	ITPSBO50	ITPSBO50W

Multiunit Scan body (2021)



Product Name WMUASB48P


Encoded Scan body (2025)

IS Encoded Scan Body		Line-Up	
		Ø4.6	Ø5.3
1	Cuff 5.0		
	型號	ISESB465S	ISESB535S
2	Cuff 7.0		
	型號	ISESB467S	ISESB537S

IT Encoded Scan Body		Line-Up	
		Ø5.8 (R/N)	Ø7.2 (W/N)
1	Cuff 2.0		
	型號	ITESB582RS	ITESB722WS
2	Cuff 4.0		
	型號	ITESB584RS	ITESB724WS

Ti-Block (Titanium Pre-Milled Blank)

Complete compatibility with the Neo Fixture.
Higher strength material incl.
screw 2 types of diameter – $\varnothing 10$ & $\varnothing 14$



Neo Activation mark

IS用 Ti-Block (鈦棒)

	Type	$\varnothing 10$	$\varnothing 14$
Neo IS	Hex	ISBKH10NBD	ISBKH14NBD
	Non-Hex	ISBKN10NBD	

Titanium Link Abutment

Titanium Grade 5 (Ti-6Al-4V)
 Titanium Link for Hybrid custom abutment Favorable to anterior tooth Abutment Screw 1ea included



IS用 Ti-Link

* Narrow / Regular / Wide / S-wide

• Diameter Ø4.0



Cuff(mm)	1.0	
Length(mm)	3.5	
Type	Hex	Non-Hex
Product name	ISLKH4035NB	ISLKN4035NB

• Diameter Ø4.5



Cuff(mm)	1.0	
Length(mm)	3.5	
Type	Hex	Non-Hex
Product name	ISLKH4535NB	ISLKN4535NB

IT用 Ti-Link



Implant	Type	Ø5.5
IT	Octa	ITLKO5535NB
	Non-Octa	ITLKN5535NB

Ti-Cylinder Abutment

Titanium Grade 5 (Ti-6Al-4V)
 Titanium Cylinder link for Hybrid custom abutment
 Abutment Screw 1ea included
 Practicable control the height by grinding
 Favorable to posterior tooth



IS用 Ti-Cylinder



Implant	Type	Ø4.5	Ø5.5
IS	Hex	ISTCAH4509	ISTCAH5509
	Non-Hex	ISTCAN4509	ISTCAN5509

Lab Analog

IS用 Lab Analog



Diameter (mm)	Ø4.1	Ø4.5
Connection	2.1 Hex	2.5 Hex
Length (mm)	8.7	8.5
Type	S-Narrow	Regular, Wide
Product name	DLAISSN	DLAISRW

* Lab Analog Screw(DLAS) included

IT用 Lab Analog




Diameter(mm)	Ø 4.8	Ø 6.5
Connection	3.1 Octa	3.1 Octa
Length(mm)	8.5	8.5
Product name	DLAITRN	DLAITWN

* Lab Analog Screw(DLAS) included


Surgical Guide Sleeve

Surgical Guide Sleeve - Narrow




Product Name	SGS55N
--------------	--------

Surgical Guide Sleeve –Regular/ Wide



Product Name	SGS70RW
--------------	---------

Surgical Guide C-Sleeve –Regular/ Wide



Product Name	CSGS70RW
--------------	----------